


Acoustic Startle and Prepulse Inhibition

Introduction and Basics

Acoustic Startle Reflex

When we hear a loud, abrupt and unexpected sound, we startle, a phenomenon that is common to all mammals [4]. The response consists of a rapid and involuntary extension and then flexion of a series of muscles. Although the startle response is a reflex, it can be modulated by many different stimuli. Startle reactions tend to be higher in the presence of threats, fear and pain [5-11] and are conversely decreased by anxiolytics [6, 12-14].




The pharmacology and circuitry of this reflex is virtually identical in rodents and humans [15-23].

Startle response is measured in rodents using automated startle chambers or “stabilimeter chambers”, with detectors recording whole-body reaction. In humans, the movements of oculomotor muscles (eye-blink reflex assessed using electromyographic recording of orbicularis oculi muscle) are most typically used.



The picture below shows still photos of the actual startle response in a rat [3].




From [3]

Acoustic Startle Circuit

In humans, changes in electrical activity in neck muscles can occur within 9 msec after the onset of an auditory stimulus and within 14 msec in jaw muscles. In rats, startle occurs within 5 msec in the neck and 8 msec in the hindleg. The very short latency of the acoustic startle reflex (ASR) indicates that a simple neural pathway must mediate it. In 1982, Davis et al [24] proposed that acoustic startle in the rat was mediated by 3 (or at most 4) synapses - 1) cochlear root neurons; 2) neurons in the nucleus reticularis pontis caudalis and 3) motoneurons in the facial motor nucleus (pinna reflex and perhaps the eyeblink) or spinal cord (whole body startle) see also [4, 11, 18, 21, 25-30]

However, this simple neural pathway may not be able to account for the relatively long latency of the eyeblink reflex (35 msec) and for this reason other pathways have been suggested.

The cochlear root neurons (in both rodents and humans) are a small group of very large cells (35 μm in diameter) embedded in the cochlear nerve. These neurons receive direct input from the spiral ganglion cells in the cochlea, making them the first acoustic neurons in the central nervous system.




Tactile Stimuli (such as a sudden air puff) can also induce a startle reaction [29] and can integrate with auditory stimuli [3, 15].

[11, 18, 21, 26-29]

Prepulse Inhibition (PPI) – Sensorimotor Gating

Prepulse Inhibition (PPI) is a neurological phenomenon in which a weak prestimulus (prepulse) inhibits the reaction to a subsequent strong, startling stimulus (pulse) [31-37]. The stimuli are usually acoustic, but tactile and visual (light and airpuff) prepulse stimuli are also used. The reduction of the startle amplitude reflects the ability of the nervous system to temporarily adapt to a strong sensory stimulus when a preceding weaker



signal is given and has also been called sensorimotor gating. PPI is detectable in numerous species ranging from mice to humans. Different prepulse-to-pulse intervals, or lead intervals, or interstimulus intervals are used ranging from 20-480 msec. Lead interval counts from the start of prepulse to the start of the pulse (onset to onset interval). Typical durations are about 20 ms for prepulse and 40 ms for pulse. Prepulse is typically set 3-16 dB louder than background. Maximum inhibition is typically observed at 40-120 ms interval in rodents, but this may vary.

Deficits of prepulse inhibition (and ASR) have been demonstrated in a wide range of neurodegenerative, psychiatric and affective disorders [23, 38-64]

Prepulse Facilitation (PPF)

The converse of prepulse inhibition, known as prepulse facilitation occurs when increased response to the target stimulus is noted after prior presentation of a non-startling stimulus. This generally occurs when the lead time (or inter-stimulus interval) between the prepulse and the target pulse exceeds 500 ms.

PPI Circuit

The circuitry and pharmacology of PPI is more complex than those of the basic ASR [2, 28, 65-69]. These are also much more controversial. Below are 2 samples of the main brain regions thought to be involved in PPI from 2 different groups.





FIGURE 10 | A schematic diagram of the startle reflex neural pathways. Adapted from [2] and [4].

Procedure

Body weight

The animals should first be weighed. Grossly different body weights can confound both ASR and PPI since the force of the muscle reaction is what is being measured. Analysis of startle by weight or inclusion of weight as a co-variant or other form of weighting/controlling for means is necessary to ensure valid results.



Acclimatization

It is important to acclimatize the subjects to the test environment [70]. The amount of time required depends on the species strain and many other variables and can only be precisely determined empirically but 5-15 min is a good estimate.

Habituation

Several (5-15) pulses of the maximum to be used (usually 110 or 120 db) should be presented [71-73]. Generally these pulses are about 40 msec in duration.


This serves several purposes. Firstly, the amount and rate of habituation (or facilitation in rare cases) are important measures. Secondly, it serves to ensure a stable baseline and it reduces variability. All sound stimuli should be presented with random time intervals (inter-trial interval, or ITI) between each presentation, or the subjects will rapidly learn to predict the impending sound. This point is somewhat controversial, as many groups provide an ITI consistent (and thus predictable)




General refs [1, 16, 17, 49, 52, 74-87]

Threshold

Next, the startle threshold should be established. A “no stim” (or 0 dB, or baseline) trial should always be included to ensure that movement, weight or other confound is not altering the apparent startle response. At least 3-5 instances of each dB levels should be presented and averaged. These should be presented in random order, with random intervals. These should NEVER be presented in ascending or descending order. Generally these pulses are about 40 msec in duration – the same values that you use during your habituation and PPI.



General refs [88-90]



PPI

There are several parameters that can be varied [91-93]. The duration and intensity (dB) of the target pulse will impact the PPI results. The inter-stimulus interval (ISI) between the prepulse and the target pulse is also a factor.

Each prepulse should be presented at least 3-6 times. Be sure to include the same number of presentations of “no stim” (0 db target pulse) and the target pulse (i.e. 110 dB) alone.

Calculating PPI

$$\%PPI = [\text{pulse-alone} - (\text{prepulse} + \text{pulse score})] / \text{pulse-alone score}$$


General refs [22, 28, 65, 67, 69, 79, 91, 94-107]

Sample PPI protocol

Pulse alone trials (P) consisted of a single white noise burst (120 dB, 40ms). The prepulse + pulse trials (such as PP69P, PP73P, PP81P) consisting of a prepulse of noise (20 ms at 69, 73, or 81 dB respectively) followed 100 ms after prepulse onset by a noise pulse (120 dB, 40 ms).

No-stimulus (NS) trials consisted of background noise only. Sessions were structured as follows: 1) 10- min acclimation at background noise level; 2) five P trials; 3) ten blocks each containing all five trials (P, PP69P, PP73P, PP81P, NS) in pseudorandom order; 4) five P trials. Inter-trials intervals were pseudorandomly distributed between 12-30 sec. The maximum force intensity for each trial was recorded as startle level. The average percent reduction in startle intensity between pulse and prepulse +pulse trials at all three prepulse levels was defined as the PPI level. The percentage PPI induced by each prepulse

intensity is calculated as $(100 - [100 \times \text{startle amplitude on pre-pulse} / \text{startle amplitude on pulse alone}])$.



Determination of target pulse.


Generally, the target pulse is around 110-120 dB, but pulses as low as 90 dB have been used. Establishing the ASR threshold first ensures that you will not have a floor effect (i.e. that startle response to the target pulse is already too low to see an effect).

Inter-stimulus Interval (ISI)

Generally an ISI of 60-200 msec has been used, and in some cases even greater ranges. Generally, an ISI of between 40-100 msec produces the maximal PPI in mice [108] (however strains are likely to vary). Be aware of the optimum range for your species, strain etc from the literature or determine it empirically. Your ability to detect differences between groups may depend on this factor, since the most robust PPI may result in a ceiling effect (or floor effect, depending on how you think about it). Longer intervals (500 msec or longer) will generally result in facilitation. [109]

Prepulse duration

The prepulse duration should be from 20-40 msec – there is no evidence that changing prepulse duration within that range alters PPI (or PPF) [110]



Reviews

[4, 111, 112]

Measuring the ASR and PPI

[113, 114]

Factors Affecting Startle and PPI

Sex

[47, 115-123]

Gonadal Hormones / Estrous/ Menstrual Cycle

[117, 120, 124-130]

Age and Development

[79, 119, 121, 131-136]

Strain

[116, 132, 137-148]

Stress and other Hormones

[30, 139, 148-150]

Ambient lighting

[41, 151-155]

Hearing Loss

[156-160]

Circadian

[161-163]

Human ASR and PPI deficits

[23, 38-64, 164-166]

Relationship to P50 (evoked potential)

[167]

Other

[168-171]

Relationship to other disorders and traits (anxiety, depression, psychiatric attention etc)

[41, 123, 150, 164, 166, 172-182]

References

1. Simons-Weidenmaier, N.S., et al., *Synaptic depression and short-term habituation are located in the sensory part of the mammalian startle pathway*. BMC Neurosci, 2006. **7**: p. 38.
2. Fendt, M., L. Li, and J.S. Yeomans, *Brain stem circuits mediating prepulse inhibition of the startle reflex*. Psychopharmacology (Berl), 2001. **156**(2-3): p. 216-24.
3. Yeomans, J.S., et al., *Tactile, acoustic and vestibular systems sum to elicit the startle reflex*. Neurosci Biobehav Rev, 2002. **26**(1): p. 1-11.
4. Koch, M., *The neurobiology of startle*. Prog Neurobiol, 1999. **59**(2): p. 107-28.
5. Davis, M., *Sensitization of the acoustic startle reflex by footshock*. Behav Neurosci, 1989. **103**(3): p. 495-503.
6. Davis, M., et al., *Fear-potentiated startle: A neural and pharmacological analysis*. Behavioural Brain Research, 1993. **58**(1-2): p. 175.
7. Hebb, A.L.O., et al., *Exposure of mice to a predator odor increases acoustic startle but does not disrupt the rewarding properties of VTA intracranial self-stimulation*. Brain Research, 2003. **982**(2): p. 195.
8. Hammond, G.R. and D.S. Leitner, *Augmentation of the rat's acoustic startle reflex by nonreflexogenic stimuli*. Behavioral Neuroscience, 1990. **104**(6): p. 841-848.
9. Kaviani, H., et al., *Modulation of the acoustic startle reflex by emotionally-toned film-clips*. International Journal of Psychophysiology, 1999. **32**(1): p. 47.
10. Sorenson, C.A. and N.R. Swerdlow, *The effect of tail pinch on the acoustic startle response in rats*. Brain Res, 1982. **247**(1): p. 105-13.
11. Davis, M., D.L. Walker, and Y. Lee, *Amygdala and bed nucleus of the stria terminalis: differential roles in fear and anxiety measured with the acoustic startle reflex*. Philos Trans R Soc Lond B Biol Sci, 1997. **352**(1362): p. 1675-87.
12. Commissaris, R.L., E.A. Fomum, and B.J. Leavell, *Effects of buspirone and alprazolam treatment on the startle-potentiated startle response*. Depress Anxiety, 2004. **19**(3): p. 146-51.
13. Riba, J., et al., *Differential effects of alprazolam on the baseline and fear-potentiated startle reflex in humans: a dose-response study*. Psychopharmacology (Berl), 2001. **157**(4): p. 358-67.
14. Rodriguez-Fornells, A., et al., *Effects of alprazolam on the acoustic startle response in humans*. Psychopharmacology (Berl), 1999. **143**(3): p. 280-5.
15. Caesar, M., J. Ostwald, and P.K. Pilz, *Startle responses measured in muscles innervated by facial and trigeminal nerves show common modulation*. Behav Neurosci, 1989. **103**(5): p. 1075-81.
16. Groves, P.M., C.J. Wilson, and R.D. Boyle, *Brain stem pathways, cortical modulation, and habituation of the acoustic startle response*. Behav Biol, 1974. **10**(4): p. 391-418.
17. Groves, P.M., C.J. Wilson, and S.W. Miller, *Habituation of the acoustic startle response: a neural systems analysis of habituation in the intact animal*. Adv Psychobiol, 1976. **3**: p. 327-79.
18. Koch, M. and H.-U. Schnitzler, *The acoustic startle response in rats--circuits mediating evocation, inhibition and potentiation*. Behavioural Brain Research, 1997. **89**(1-2): p. 35.
19. Kumari, V., et al., *Structural brain correlates of prepulse inhibition of the acoustic startle response in healthy humans*. NeuroImage, 2005. **26**(4): p. 1052.
20. Reijmers, L.G.J.E. and B.W.M.M. Peeters, *Acoustic prepulses can facilitate the startle reflex in rats: Discrepancy between rat and human data resolved*. Brain Research Bulletin, 1994. **35**(4): p. 337-338.
21. Yeomans, J.S. and P.W. Frankland, *The acoustic startle reflex: neurons and connections*. Brain Research Reviews, 1995. **21**(3): p. 301.
22. Chien, A.S., et al., *Measuring prepulse inhibition and facilitation of the eyeblink response to acoustic startle during continuous acquisition with ER-fMRI*. Schizophrenia Research, 2003. **60**(1, Supplement 1): p. 214-215.

23. Braff, D., M. Geyer, and N. Swerdlow, *Human studies of prepulse inhibition of startle: normal subjects, patient groups, and pharmacological studies*. Psychopharmacology, 2001. **156**(2): p. 234-258.
24. Davis, M., et al., *A primary acoustic startle circuit: lesion and stimulation studies*. J Neurosci, 1982. **2**(6): p. 791-805.
25. Lee, Y., et al., *A primary acoustic startle pathway: obligatory role of cochlear root neurons and the nucleus reticularis pontis caudalis*. J Neurosci, 1996. **16**(11): p. 3775-89.
26. Davis, M., *Pharmacological and anatomical analysis of fear conditioning using the fear-potentiated startle paradigm*. Behavioral Neuroscience, 1986. **100**(6): p. 814-824.
27. Plappert, C.F. and P.K. Pilz, *The acoustic startle response as an effective model for elucidating the effect of genes on the neural mechanism of behavior in mice*. Behav Brain Res, 2001. **125**(1-2): p. 183-8.
28. Swerdlow, et al., *Neural circuit regulation of prepulse inhibition of startle in the rat: current knowledge and future challenges*. Psychopharmacology, 2001. **156**(2): p. 194-215.
29. Schmid, S., N.S. Simons, and H.U. Schnitzler, *Cellular mechanisms of the trigeminally evoked startle response*. Eur J Neurosci, 2003. **17**(7): p. 1438-44.
30. Lee, Y. and M. Davis, *Role of the hippocampus, the bed nucleus of the stria terminalis, and the amygdala in the excitatory effect of corticotropin-releasing hormone on the acoustic startle reflex*. J Neurosci, 1997. **17**(16): p. 6434-46.
31. Hoffman, H.S. and M. Fleshler, *Startle Reaction: Modification By Background Acoustic Stimulation*. Science, 1963. **141**: p. 928-30.
32. Hoffman, H.S., R.R. Marsh, and N. Stein, *Persistence of background acoustic stimulation in controlling startle*. J Comp Physiol Psychol, 1969. **68**(2): p. 280-3.
33. Hoffman, H.S. and J.L. Searle, *Acoustic Variables In The Modification Of Startle Reaction In The Rat*. J Comp Physiol Psychol, 1965. **60**: p. 53-8.
34. Hoffman, H.S. and B.L. Wible, *Role of weak signals in acoustic startle*. J Acoust Soc Am, 1970. **47**(2): p. 489-97.
35. Ison, J.R., D.W. McAdam, and G.R. Hammond, *Latency and amplitude changes in the acoustic startle reflex of the rat produced by variation in auditory prestimulation*. Physiol Behav, 1973. **10**(6): p. 1035-9.
36. Ison, J.R., et al., *Facilitation and Inhibition of the Acoustic Startle Reflex in the Rat After a Momentary Increase in Background Noise Level*. Behavioral Neuroscience, 1997. **111**(6): p. 1335.
37. Marsh, R.R., et al., *The role of small changes in the acoustic environment in modifying the startle reflex*. J Exp Psychol Anim Behav Process, 1975. **1**(3): p. 235-44.
38. Ford, J.M., et al., *Automatic and Effortful Processing in Aging and Dementia: Event-Related Brain Potentials*. Neurobiology of Aging, 1997. **18**(2): p. 169-180.
39. Geyer, M.A., et al., *Isolation rearing of rats produces a deficit in prepulse inhibition of acoustic startle similar to that in schizophrenia*. Biol Psychiatry, 1993. **34**(6): p. 361-72.
40. Gironell, A., et al., *Abnormalities of the acoustic startle reflex and reaction time in gilles de la tourette syndrome*. Clin Neurophysiol, 2000. **111**(8): p. 1366-71.
41. Grillon, C., et al., *Effect of darkness on acoustic startle in Vietnam veterans with PTSD*. Am J Psychiatry, 1998. **155**(6): p. 812-7.
42. Guitar, B., *Acoustic startle responses and temperament in individuals who stutter*. J Speech Lang Hear Res, 2003. **46**(1): p. 233-40.
43. Hessel, D., et al., *Amygdala dysfunction in men with the fragile X premutation*. Brain, 2007. **130**(Pt 2): p. 404-16.
44. Hoenig, K., et al., *Impaired Prepulse Inhibition of Acoustic Startle in Obsessive-Compulsive Disorder*. Biological Psychiatry, 2005. **57**(10): p. 1153.
45. Jankelowitz, S.K. and J.G. Colebatch, *The acoustic startle reflex in ischemic stroke*. Neurology, 2004. **62**(1): p. 114-6.
46. Klorman, R., et al., *Acoustic startle in maltreated children*. J Abnorm Child Psychol, 2003. **31**(4): p. 359-70.
47. Kumari, V., I. Aasen, and T. Sharma, *Sex differences in prepulse inhibition deficits in chronic schizophrenia*. Schizophrenia Research, 2004. **69**(2-3): p. 219-235.

48. Lipschitz, D.S., et al., *Baseline and modulated acoustic startle responses in adolescent girls with posttraumatic stress disorder*. J Am Acad Child Adolesc Psychiatry, 2005. **44**(8): p. 807-14.
49. Ludewig, K., et al., *Stability of the acoustic startle reflex, prepulse inhibition, and habituation in schizophrenia*. Schizophrenia Research, 2002. **55**(1-2): p. 129.
50. Ludewig, K., M.A. Geyer, and F.X. Vollenweider, *Deficits in prepulse inhibition and habituation in never-medicated, first-episode schizophrenia*. Biological Psychiatry, 2003. **54**(2): p. 121-128.
51. Ludewig, S., et al., *Information-processing deficits and cognitive dysfunction in panic disorder*. J Psychiatry Neurosci, 2005. **30**(1): p. 37-43.
52. Meincke, U., et al., *Sensitization and habituation of the acoustic startle reflex in patients with schizophrenia*. Psychiatry Research, 2004. **126**(1): p. 51.
53. Morgan, C.A., 3rd, et al., *Exaggerated acoustic startle reflex in Gulf War veterans with posttraumatic stress disorder*. Am J Psychiatry, 1996. **153**(1): p. 64-8.
54. Parwani, A., et al., *Impaired prepulse inhibition of acoustic startle in schizophrenia*. Biological Psychiatry, 2000. **47**(7): p. 662.
55. Sanchez-Navarro, J.P., J.M. Martínez-Selva, and F. Román, *Emotional Response in Patients With Frontal Brain Damage: Effects of Affective Valence and Information Content*. Behavioral Neuroscience, 2005. **119**(1): p. 87-97.
56. Siegelar, S.E., et al., *The auditory startle response in post-traumatic stress disorder*. Exp Brain Res, 2006. **174**(1): p. 1-6.
57. Swerdlow, N.R., et al., *Impaired prepulse inhibition of acoustic and tactile startle response in patients with Huntington's disease*. J Neurol Neurosurg Psychiatry, 1995. **58**(2): p. 192-200.
58. Wynn, J.K., et al., *Prepulse facilitation deficits in a passive task in schizophrenia patients and their siblings*. Schizophrenia Research, 2003. **60**(1, Supplement 1): p. 264-215.
59. Wynn, J.K., et al., *Prepulse facilitation and prepulse inhibition in schizophrenia patients and their unaffected siblings*. Biological Psychiatry, 2004. **55**(5): p. 518-523.
60. Wynn, J.K., et al., *Sensorimotor gating, orienting and social perception in schizophrenia*. Schizophrenia Research, 2005. **73**(2-3): p. 319-325.
61. Forbes, E.E., et al., *Affect-modulated startle in adults with childhood-onset depression: relations to bipolar course and number of lifetime depressive episodes*. Psychiatry Res, 2005. **134**(1): p. 11-25.
62. Kaviani, H., et al., *Affective modulation of the startle response in depression: influence of the severity of depression, anhedonia, and anxiety*. J Affect Disord, 2004. **83**(1): p. 21-31.
63. Perry, W., A. Minassian, and D. Feifel, *Prepulse inhibition in patients with non-psychotic major depressive disorder*. J Affect Disord, 2004. **81**(2): p. 179-84.
64. Ludewig, S. and K. Ludewig, *No prepulse inhibition deficits in patients with unipolar depression*. Depress Anxiety, 2003. **17**(4): p. 224-5.
65. Caine, S.B., M.A. Geyer, and N.R. Swerdlow, *Hippocampal modulation of acoustic startle and prepulse inhibition in the rat*. Pharmacology Biochemistry and Behavior, 1992. **43**(4): p. 1201.
66. Fendt, M., M. Koch, and H.-U. Schnitzler, *Lesions of the central gray block the sensitization of the acoustic startle response in rats*. Brain Research, 1994. **661**(1-2): p. 163.
67. Goto, K., et al., *Reduced prepulse inhibition in rats with entorhinal cortex lesions*. Behav Brain Res, 2002. **134**(1-2): p. 201-7.
68. Goto, K., et al., *Involvement of nucleus accumbens dopaminergic transmission in acoustic startle: observations concerning prepulse inhibition in rats with entorhinal cortex lesions*. Psychiatry Clin Neurosci, 2004. **58**(4): p. 441-5.
69. Koch, M., *The septohippocampal system is involved in prepulse inhibition of the acoustic startle response in rats*. Behav Neurosci, 1996. **110**(3): p. 468-77.
70. Faraday, M.M. and N.E. Grunberg, *The Importance of Acclimation in Acoustic Startle Amplitude and Pre-Pulse Inhibition Testing of Male and Female Rats*. Pharmacology Biochemistry and Behavior, 2000. **66**(2): p. 375.

71. Pilz, P.K.D. and H.-U. Schnitzler, *Habituation and Sensitization of the Acoustic Startle Response in Rats: Amplitude, Threshold, and Latency Measures*. Neurobiology of Learning and Memory, 1996. **66**(1): p. 67.
72. Pilz, P.K., H.U. Schnitzler, and D. Menne, *Acoustic startle threshold of the albino rat (*Rattus norvegicus*)*. J Comp Psychol, 1987. **101**(1): p. 67-72.
73. Pinckney, L. and J.R. Ison, *Modification of primary and secondary startle reactions to acoustic stimuli in the rat*. Psychophysiology, 1979. **16**(5): p. 420-8.
74. Storozeva, Z.I. and M.V. Pletnicov, *Habituation of acoustic startle in rats--a functional ablation study*. Neuroreport, 1994. **5**(16): p. 2065-8.
75. Weber, M., H.U. Schnitzler, and S. Schmid, *Synaptic plasticity in the acoustic startle pathway: the neuronal basis for short-term habituation?* Eur J Neurosci, 2002. **16**(7): p. 1325-32.
76. Wilson, C.J. and P.M. Groves, *Refractory period and habituation of acoustic startle response in rats*. J Comp Physiol Psychol, 1973. **83**(3): p. 492-8.
77. Bespalov, A., et al., *Habituation deficits induced by metabotropic glutamate receptors 2/3 receptor blockade in mice: reversal by antipsychotic drugs*. J Pharmacol Exp Ther, 2007. **320**(2): p. 944-50.
78. Cassella, J.V. and M. Davis, *Habituation, prepulse inhibition, fear conditioning, and drug modulation of the acoustically elicited pinna reflex in rats*. Behavioral Neuroscience, 1986. **100**(1): p. 39-44.
79. Ellwanger, J., M.A. Geyer, and D.L. Braff, *The relationship of age to prepulse inhibition and habituation of the acoustic startle response*. Biological Psychology, 2003. **62**(3): p. 175.
80. Hoffman, H.S. and C. Stitt, *Behavioral factors in habituation of acoustic startle reactions*. J Comp Physiol Psychol, 1969. **68**(2): p. 276-9.
81. Leaton, R.N. and W.F. Supple, *Medial cerebellum and long-term habituation of acoustic startle in rats*. Behavioral Neuroscience, 1991. **105**(6): p. 804-816.
82. Pilz, P.K., T.D. Carl, and C.F. Plappert, *Habituation of the acoustic and the tactile startle responses in mice: two independent sensory processes*. Behav Neurosci, 2004. **118**(5): p. 975-83.
83. Pilz, P.K. and H.U. Schnitzler, *Habituation and sensitization of the acoustic startle response in rats: amplitude, threshold, and latency measures*. Neurobiol Learn Mem, 1996. **66**(1): p. 67-79.
84. Pilz, P.K.D., T.D. Carl, and C.F. Plappert, *Habituation of the Acoustic and the Tactile Startle Responses in Mice: Two Independent Sensory Processes*. Behavioral Neuroscience, 2004. **118**(5): p. 975-983.
85. Plappert, C.F. and P.K. Pilz, *Long-term habituation of the startle response in mice evoked by acoustic and tactile stimuli*. Behav Brain Res, 2005. **162**(2): p. 307-10.
86. Sanford, L.D., et al., *Peripheral and Central Components of Alerting: Habituation of Acoustic Startle, Orienting Responses, and Elicited Waveforms*. Behavioral Neuroscience, 1992. **106**(1): p. 112.
87. Davis, M. and A.R. Wagner, *Habituation of startle response under incremental sequence of stimulus intensities*. J Comp Physiol Psychol, 1969. **67**(4): p. 486-92.
88. Blumenthal, T.D., *The startle response to acoustic stimuli near startle threshold: effects of stimulus rise and fall time, duration, and intensity*. Psychophysiology, 1988. **25**(5): p. 607-11.
89. Pilz, P.K., M. Caeser, and J. Ostwald, *Comparative threshold studies of the acoustic pinna, jaw and startle reflex in the rat*. Physiol Behav, 1988. **43**(4): p. 411-5.
90. Davis, M., *Signal-to-noise ratio as a predictor of startle amplitude and habituation in the rat*. J Comp Physiol Psychol, 1974. **86**(5): p. 812-25.
91. Plappert, C.F., P.K. Pilz, and H.U. Schnitzler, *Factors governing prepulse inhibition and prepulse facilitation of the acoustic startle response in mice*. Behav Brain Res, 2004. **152**(2): p. 403-12.
92. Reijmers, L.G. and B.W. Peeters, *Effects of acoustic prepulses on the startle reflex in rats: a parametric analysis*. Brain Res, 1994. **661**(1-2): p. 174-80.

93. Sasaki, H., et al., *Prepulse facilitation of auditory startle response in hamsters*. Neuroscience Letters, 1998. **248**(2): p. 117-120.
94. Yee, B.K., et al., *The expression of prepulse inhibition of the acoustic startle reflex as a function of three pulse stimulus intensities, three prepulse stimulus intensities, and three levels of startle responsiveness in C57BL6/J mice*. Behavioural Brain Research, 2005. **163**(2): p. 265.
95. Bitsios, P. and S.G. Giakoumaki, *Relationship of prepulse inhibition of the startle reflex to attentional and executive mechanisms in man*. Int J Psychophysiol, 2005. **55**(2): p. 229-41.
96. Blumenthal, T.D. and G.A. Gescheider, *Modification of the acoustic startle reflex by a tactile prepulse: the effects of stimulus onset asynchrony and prepulse intensity*. Psychophysiology, 1987. **24**(3): p. 320-7.
97. Csomor, P.A., et al., *On the feasibility to detect and to quantify prepulse-elicited reaction in prepulse inhibition of the acoustic startle reflex in humans*. Behavioural Brain Research, 2005. **162**(2): p. 256.
98. Csomor, P.A., et al., *The monotonic dependency of prepulse inhibition of the acoustic startle reflex on the intensity of the startle-eliciting stimulus*. Behavioural Brain Research, 2006. **174**(1): p. 143-150.
99. Gewirtz, J.C. and M. Davis, *Habituation of prepulse inhibition of the startle reflex using an auditory prepulse close to background noise*. Behavioral Neuroscience, 1995. **109**(3): p. 388-395.
100. Johansson, C., et al., *Prepulse inhibition of acoustic startle, a measure of sensorimotor gating: Effects of antipsychotics and other agents in rats*. Pharmacology Biochemistry and Behavior, 1995. **52**(4): p. 649.
101. Kodsi, M.H. and N.R. Swerdlow, *Ventral pallidal GABA-A receptors regulate prepulse inhibition of acoustic startle*. Brain Res, 1995. **684**(1): p. 26-35.
102. Kodsi, M.H. and N.R. Swerdlow, *Prepulse inhibition in the rat is regulated by ventral and caudodorsal striato-pallidal circuitry*. Behavioral Neuroscience, 1995. **109**(5): p. 912-928.
103. Leumann, L., et al., *A neural network approach to the acoustic startle reflex and prepulse inhibition*. Brain Res Bull, 2001. **56**(2): p. 101-10.
104. Meincke, U., et al., *On the waveform of the acoustic startle blink in the paradigm of prepulse inhibition - methodological and physiological aspects*. Neuropsychobiology, 2005. **52**(1): p. 24-7.
105. Pryce, C.R., et al., *Comparison of the effects of infant handling, isolation, and nonhandling on acoustic startle, prepulse inhibition, locomotion, and HPA activity in the adult rat*. Behav Neurosci, 2001. **115**(1): p. 71-83.
106. Schwarzkopf, S.B., et al., *Test-retest reliability of prepulse inhibition of the acoustic startle response*. Biological Psychiatry, 1993. **34**(12): p. 896.
107. Takahashi, K., et al., *Neural Circuits Containing Pallidotegmental GABAergic Neurons are Involved in the Prepulse Inhibition of the Startle Reflex in Mice*. Biological Psychiatry, 2007. **62**(2): p. 148-157.
108. Plappert, C.F., P.K.D. Pilz, and H.-U. Schnitzler, *Factors governing prepulse inhibition and prepulse facilitation of the acoustic startle response in mice*. Behavioural Brain Research, 2004. **152**(2): p. 403.
109. Lipp, O.V. and D.A.T. Siddle, *The effects of prepulse-blink reflex trial repetition and prepulse change on blink reflex modification at short and long lead intervals*. Biological Psychology, 1998. **47**(1): p. 45-63.
110. Reijmers, L.G.J.E. and B.W.M.M. Peeters, *Effects of acoustic prepulses on the startle reflex in rats: a parametric analysis*. Brain Research, 1994. **667**(1): p. 144-150.
111. Grillon, C. and J. Baas, *A review of the modulation of the startle reflex by affective states and its application in psychiatry*. Clin Neurophysiol, 2003. **114**(9): p. 1557-79.
112. Grillon, C., *Startle reactivity and anxiety disorders: aversive conditioning, context, and neurobiology*. Biol Psychiatry, 2002. **52**(10): p. 958-75.
113. Light, G.A., et al., *Stability of the acoustic startle reflex: Amplitude vs area measures*. Biological Psychiatry, 1995. **37**(9): p. 618.

114. Walker, D.L. and M. Davis, *Quantifying fear potentiated startle using absolute versus proportional increase scoring methods: implications for the neurocircuitry of fear and anxiety*. Psychopharmacology (Berl), 2002. **164**(3): p. 318-28.
115. Aasen, I., L. Kolli, and V. Kumari, *Sex effects in prepulse inhibition and facilitation of the acoustic startle response: implications for pharmacological and treatment studies*. J Psychopharmacol, 2005. **19**(1): p. 39-45.
116. Blaszczyk, J. and K. Tajchert, *Sex and strain differences of acoustic startle reaction development in adolescent albino Wistar and hooded rats*. Acta Neurobiol Exp (Wars), 1996. **56**(4): p. 919-25.
117. Ison, J.R. and P.D. Allen, *Pre- but not post-menopausal female CBA/CaJ mice show less prepulse inhibition than male mice of the same age*. Behavioural Brain Research, 2007. **185**(2): p. 76-81.
118. Kumari, V., et al., *Sex differences in prepulse inhibition of the acoustic startle response*. Personality and Individual Differences, 2003. **35**(4): p. 733.
119. Ludewig, K., et al., *The acoustic startle reflex and its modulation: effects of age and gender in humans*. Biological Psychology, 2003. **63**(3): p. 311.
120. Plappert, C.F., A.M. Rodenbucher, and P.K.D. Pilz, *Effects of sex and estrous cycle on modulation of the acoustic startle response in mice*. Physiology & Behavior, 2005. **84**(4): p. 585-594.
121. van den Buuse, M., E.R. Simpson, and M.E. Jones, *Prepulse inhibition of acoustic startle in aromatase knock-out mice: effects of age and gender*. Genes Brain Behav, 2003. **2**(2): p. 93-102.
122. Toufexis, D., et al., *Sex differences in hormonal modulation of anxiety measured with light-enhanced startle: possible role for arginine vasopressin in the male*. J Neurosci, 2005. **25**(39): p. 9010-6.
123. Adamec, R., et al., *Lasting anxiogenic effects of feline predator stress in mice: sex differences in vulnerability to stress and predicting severity of anxiogenic response from the stress experience*. Physiol Behav, 2006. **88**(1-2): p. 12-29.
124. Jovanovic, T., et al., *Menstrual cycle phase effects on prepulse inhibition of acoustic startle*. Psychophysiology, 2004. **41**(3): p. 401-6.
125. Karper, L.P., et al., *Pre-pulse inhibition of the acoustic startle response and measures of attention in schizophrenic subjects*. Schizophrenia Research, 1993. **9**(2-3): p. 160.
126. Swerdlow, N.R., P.L. Hartman, and P.P. Auerbach, *Changes in sensorimotor inhibition across the menstrual cycle: Implications for neuropsychiatric disorders*. Biological Psychiatry, 1997. **41**(4): p. 452-460.
127. Turvin, J.C., J.W.S. Messer, and M.F. Kritzer, *On again, off again effects of gonadectomy on the acoustic startle reflex in adult male rats*. Physiology & Behavior, 2007. **90**(2-3): p. 473-482.
128. Vaillancourt, C., et al., *Effects of ovariectomy and estradiol on acoustic startle responses in rats*. Pharmacology Biochemistry and Behavior, 2002. **74**(1): p. 103.
129. Van den Buuse, M. and N. Eikelis, *Estrogen increases prepulse inhibition of acoustic startle in rats*. European Journal of Pharmacology, 2001. **425**(1): p. 33.
130. Koch, M., *Sensorimotor gating changes across the estrous cycle in female rats*. Physiol Behav, 1998. **64**(5): p. 625-8.
131. Acri, J.B., et al., *Strain and age differences in acoustic startle responses and effects of nicotine in rats*. Pharmacology Biochemistry and Behavior, 1995. **50**(2): p. 191.
132. Acri, J.B., et al., *Effects of age, strain, and nicotine on rats' acoustic startle*. Pharmacology Biochemistry and Behavior, 1993. **46**(2): p. 505.
133. Barsz, K., et al., *Behavioral and neural measures of auditory temporal acuity in aging humans and mice*. Neurobiology of Aging, 2002. **23**(4): p. 565-578.
134. Cripps, M.M. and D.J. Nash, *Ontogeny and adult behavior of mice with congenital neural tube defects*. Behavioral and Neural Biology, 1983. **38**(1): p. 127-132.
135. Parham, K. and J.F. Willott, *Acoustic startle response in young and aging C57BL/6J and CBA/J mice*. Behav Neurosci, 1988. **102**(6): p. 881-6.

136. Pletnicov, M.V., Z.I. Storozheva, and V.V. Sherstnev, *Developmental analysis of habituation of the acoustic startle response in the preweanling and adult rats*. Behavioural Processes, 1995. **34**(3): p. 269.
137. Bullock, A.E., et al., *Inbred mouse strains differ in the regulation of startle and prepulse inhibition of the startle response*. Behavioral Neuroscience, 1997. **111**(6): p. 1353-1360.
138. Commissaris, R.L., et al., *MR/Har and MNRA/Har Maudsley rat strains: differences in acoustic startle habituation*. Behav Genet, 1988. **18**(6): p. 663-9.
139. Conti, L.H. and M.P. Printz, *Rat strain-dependent effects of repeated stress on the acoustic startle response*. Behavioural Brain Research, 2003. **144**(1-2): p. 11.
140. Falls, W.A., et al., *Fear-Potentiated Startle in Two Strains of Inbred Mice*. Behavioral Neuroscience, 1997. **111**(4): p. 855-861.
141. Lehmann, J., C.R. Pryce, and J. Feldon, *Sex differences in the acoustic startle response and prepulse inhibition in Wistar rats*. Behavioural Brain Research, 1999. **104**(1-2): p. 113.
142. Paylor, R., D. Pickar, and J.N. Crawley, *Characterization of inbred mouse strain differences in startle and the prepulse inhibition model of schizophrenia*. Schizophrenia Research, 1997. **24**(1-2): p. 59.
143. Plappert, C.F. and P.K. Pilz, *Difference in anxiety and sensitization of the acoustic startle response between the two inbred mouse strains BALB/cAN and DBA/2N*. Genes Brain Behav, 2002. **1**(3): p. 178-86.
144. Swerdlow, N.R., et al., *Sensitivity to drug effects on prepulse inhibition in inbred and outbred rat strains*. Pharmacology Biochemistry and Behavior, 2004. **77**(2): p. 291.
145. van den Buuse, M., *Deficient prepulse inhibition of acoustic startle in Hooded-Wistar rats compared with Sprague-Dawley rats*. Clin Exp Pharmacol Physiol, 2003. **30**(4): p. 254-61.
146. Willott, J.F., et al., *Acoustic Startle and Prepulse Inhibition in 40 Inbred Strains of Mice*. Behavioral Neuroscience, 2003. **117**(4): p. 716.
147. Wolf, R., et al., *Prepulse inhibition is different in two inbred mouse strains (CPB-K and BALB/cJ) with different hippocampal NMDA receptor densities*. Behavioural Brain Research, 2006. **166**(1): p. 78-84.
148. Conti, L.H. and M.P. Printz, *Rat strain-dependent effects of repeated stress on the acoustic startle response*. Behav Brain Res, 2003. **144**(1-2): p. 11-8.
149. Acri, J.B., et al., *Interactions of stress and nicotine on PPI and acoustic startle*. Pharmacology Biochemistry and Behavior, 1993. **46**(2): p. 505.
150. Grillon, C., et al., *Cortisol and DHEA-S are associated with startle potentiation during aversive conditioning in humans*. Psychopharmacology (Berl), 2005: p. 1-8.
151. Basso, M.R., Jr., *Neurobiological relationships between ambient lighting and the startle response to acoustic stress in humans*. Int J Neurosci, 2001. **110**(3-4): p. 147-57.
152. Grillon, C., M. Pellowski, and K. Merikangas, *The startle reflex is facilitated by darkness in humans*. Biological Psychiatry, 1996. **39**(7): p. 555-460.
153. Grillon, C., et al., *Darkness facilitates the acoustic startle reflex in humans*. Biological Psychiatry, 1997. **42**(6): p. 453.
154. Walker, D.L. and M. Davis, *Anxiogenic effects of high illumination levels assessed with the acoustic startle response in rats*. Biological Psychiatry, 1997. **42**(6): p. 461.
155. Weber, M., N. Watts, and R. Richardson, *High illumination levels potentiate the acoustic startle response in preweanling rats*. Behav Neurosci, 2003. **117**(6): p. 1458-62.
156. Carlson, S. and J.F. Willott, *The behavioral salience of tones as indicated by prepulse inhibition of the startle response: relationship to hearing loss and central neural plasticity in C57BL/6J mice*. Hearing Research, 1996. **99**(1-2): p. 168-175.
157. Ouagazzal, A.-M., D. Reiss, and R. Romand, *Effects of age-related hearing loss on startle reflex and prepulse inhibition in mice on pure and mixed C57BL and 129 genetic background*. Behavioural Brain Research, 2006. **172**(2): p. 307-315.
158. Willott, J.F. and S. Carlson, *Modification of the Acoustic Startle Response in Hearing-Impaired C57BL/6J Mice: Prepulse Augmentation and Prolongation of Prepulse Inhibition*. Behavioral Neuroscience, 1995. **109**(3): p. 396.

159. Willott, J.F., S. Carlson, and H. Chen, *Prepulse Inhibition of the Startle Response in Mice: Relationship to Hearing Loss and Auditory System Plasticity*. Behavioral Neuroscience, 1994. **108**(4): p. 703-713.
160. Willott, J.F., et al., *The BALB/c mouse as an animal model for progressive sensorineural hearing loss*. Hearing Research, 1998. **115**(1-2): p. 162-174.
161. Chabot, C.C. and D.H. Taylor, *Circadian Modulation of the Rat Acoustic Startle Response*. Behavioral Neuroscience, 1992. **106**(5): p. 846.
162. Chabot, C.C. and D.H. Taylor, *Daily rhythmicity of the rat acoustic startle response*. Physiology & Behavior, 1992. **51**(4): p. 885.
163. Frankland, P.W. and M.R. Ralph, *Circadian modulation in the rat acoustic startle circuit*. Behav Neurosci, 1995. **109**(1): p. 43-8.
164. Hazlett, E.A., et al., *Exaggerated affect-modulated startle during unpleasant stimuli in borderline personality disorder*. Biol Psychiatry, 2007. **62**(3): p. 250-5.
165. Bernier, R., et al., *Individuals with autism spectrum disorder show normal responses to a fear potential startle paradigm*. J Autism Dev Disord, 2005. **35**(5): p. 575-83.
166. Grillon, C., L. Dierker, and K.R. Merikangas, *Startle modulation in children at risk for anxiety disorders and/or alcoholism*. J Am Acad Child Adolesc Psychiatry, 1997. **36**(7): p. 925-32.
167. Schwarzkopf, S.B., J.S. Lamberti, and D.A. Smith, *Concurrent assessment of acoustic startle and auditory P50 evoked potential measures of sensory inhibition*. Biol Psychiatry, 1993. **33**(11-12): p. 815-28.
168. Armus, H.L. and K.A. Isham, *Effect of hunger on the acoustic startle reflex*. Psychol Rep, 1984. **54**(2): p. 667-70.
169. Armus, H.L. and K.A. Isham, *Effect of thirst on the acoustic startle reflex*. Psychol Rep, 1985. **56**(1): p. 88-90.
170. Schmid, A., M. Koch, and H.U. Schnitzler, *Conditioned pleasure attenuates the startle response in rats*. Neurobiol Learn Mem, 1995. **64**(1): p. 1-3.
171. Paschall, G.Y. and M. Davis, *Olfactory-mediated fear-potentiated startle*. Behav Neurosci, 2002. **116**(1): p. 4-12.
172. Dieckmann, M., et al., *Disturbed social behavior and motivation in rats selectively bred for deficient sensorimotor gating*. Schizophrenia Research. **In Press, Corrected Proof**.
173. Freudenberg, F., et al., *Selective breeding for deficient sensorimotor gating is accompanied by increased perseveration in rats*. Neuroscience, 2007. **148**(3): p. 612-622.
174. Ardayfio, P. and K.S. Kim, *Anxiogenic-like effect of chronic corticosterone in the light-dark emergence task in mice*. Behav Neurosci, 2006. **120**(2): p. 249-56.
175. Cohen, H., et al., *Decreased circulatory levels of neuroactive steroids in behaviourally more extremely affected rats subsequent to exposure to a potentially traumatic experience*. Int J Neuropsychopharmacol, 2007. **10**(2): p. 203-9.
176. Grillon, C., J. Levenson, and D.S. Pine, *A single dose of the selective serotonin reuptake inhibitor citalopram exacerbates anxiety in humans: a fear-potentiated startle study*. Neuropsychopharmacology, 2007. **32**(1): p. 225-31.
177. Arikian, K., et al., *EEG correlates of startle reflex with reactivity to eye opening in psychiatric disorders: preliminary results*. Clin EEG Neurosci, 2006. **37**(3): p. 230-4.
178. Cornwell, B.R., et al., *Anticipation of public speaking in virtual reality reveals a relationship between trait social anxiety and startle reactivity*. Biol Psychiatry, 2006. **59**(7): p. 664-6.
179. Prehn, A., et al., *Chemosensory anxiety signals augment the startle reflex in humans*. Neurosci Lett, 2006. **394**(2): p. 127-30.
180. Harris, A.C. and J.C. Gewirtz, *Elevated startle during withdrawal from acute morphine: a model of opiate withdrawal and anxiety*. Psychopharmacology (Berl), 2004. **171**(2): p. 140-7.
181. Lissek, S. and A.S. Powers, *Sensation seeking and startle modulation by physically threatening images*. Biol Psychol, 2003. **63**(2): p. 179-97.
182. Kumari, V., et al., *Enhanced startle reactions to acoustic stimuli in patients with obsessive-compulsive disorder*. Am J Psychiatry, 2001. **158**(1): p. 134-6.

