

High throughput Forced Swim Test Analysis

Europe:
BIOBSERVE GmbH
Koenigswinterer Str. 418
53227 Bonn
Tel.: +49 228 20 160 20
Fax: +49 228 20 160 29

www.biobserve.com info@biobserve.com

North America:
BIOBSERVE
2125 Center Ave., Suite 500
Fort Lee, New Jersey 07024
Tel.: 201-302-6083
Fax: 201-302-6062



Introduction

FST Analysis



This software is designed to automate the forced swim test. The software analyzes the video image of up to eight animals simultaneously and automatically calculates the time spent struggling, swimming and floating as suggested by newly enhanced approaches.



Features I

FST Analysis

Mice & Rats

Mice and rats show different behaviors during the forced swim test. Our system can be used with both species.



High throughput

Up to 8 animals can be tested simultaneously what allows high throughput testing. All acquired data are calculated and presented in real time already during the experiment. The possibility to save different setup configurations and to generate pdf based summary reports automatically makes it even easier to deal with your test sessions.



Three behaviors

Comparing the active versus the inactive time provides only a small amount of results that are possible to get from the forced swim test. Our system is capable to distinguish between the three behaviors struggling, swimming and floating and considers also the temporal course of the shown behaviors.

Struggling
Swimming
Floating




Features II

FST Analysis

Real time individual and group performance calculation

During the experiment the acquired data are calculated and displayed graphically in real time. Thus you are able to really see what is going on.




Extensive statistic functions

The software offers a big number of statistical test to analyze the acquired data. You don't need to go through the data with other programs. Just select the algebraic function or test you want to be calculated and see the results instantly. One behavioral test - one software.

$$\begin{aligned} 2\theta_{ij} \frac{\partial P}{\partial T} &= 2\theta_{ij} \frac{P_{m,n}^{t+1} - P_{m,n}^t}{\Delta t} \\ k_s \frac{1}{r} \frac{\partial P^2}{\partial r} &= k_s \frac{P_{m+1,n}^2 - P_{m-1,n}^2}{2r\Delta r} \\ k_s \frac{\partial^2 P^2}{\partial r^2} &= k_s \frac{P_{m+1,n}^2 + P_{m-1,n}^2 - 2P_{m,n}^2}{\Delta r^2} \\ k_t \frac{\partial^2 P^2}{\partial t^2} &= k_t \frac{P_{m,n}^2 + P_{m,n}^{t+1} - 2P_{m,n}^t}{\Delta t^2} \end{aligned}$$

Useful diagrams

Already during the experiment the acquired data are analyzed and presented in convincing diagrams. All diagrams can be exported in different file formats (bmp, wmf and emf).





Features III

FST Analysis

Report generator

An integrated Acrobat pdf based report generator makes it easy to generate standardized experiment reports for each animal or group in a very short time.


Further processing of data and diagrams

The collected data can be exported for other applications. Further individual data manipulations can be conducted.

All diagrams that have been generated during the experiment or the data analysis are saved in the vector graphic format (wmf) thereby making them easy to integrate and work with in other programs.

Database integration

The software offers standard interfaces to communicate with other applications like databases. Thereby it is for example possible to transfer data into existing animal databases easily. You can import database records into the FST software, perform experiments and export the data back into your animal database. That is true high throughput!



Export filters for the data:
*.txt, *.doc, *.csv, *.xls, *.xml

Export filters for the
diagrams:
*.bmp, *.wmf, *.emf






Features IV

FST Analysis

FireWire support

The software supports the FireWire/iLink port for digital video cameras. You can use any DV-Camera with a FireWire port. This simplifies the installation because no framegrabber card or driver problems occur. It's just plug and measure.



Digital video encoding included

You do not want to fill another locker with your experiment tapes? With the FST Analysis software you can directly encode your experiment to a digital video file in real time. Easy to do, easy to store and easy to handle.



Phenotyping with the FST system

The following list shows what kind of phenotypes can be detected with the FST system based on the Mouse Genome Informatics (MGI) project that provides integrated access to data on the genetics, genomics, and biology of the laboratory mouse.

Name of the test (MP term)	MP ID
Behavioral despair	2573
decreased coping response	3064
increased coping response	3063



Program Tour - Configuration

FST Analysis

Experiment description

Experiment | Technical preferences

Description

Supervisor: Experiment ID:

Short desc.:

Free 1: Free 2:

Free 3: Free 4:

Long description:

Timing

Experiment length [s]:

Interval length [s]:

Coloring and behavior names

Struggle Set color...

Swim Set color...

Float Set color...

Calculation interval

Interval length [ms]:

Number of animals

☐ One animal
☐ Two animals
☐ Four animals
☒ Eight animals

	1	2	3	4
1 ID	1	2		
SD	Dopa 100	Dopa 200		
Group	1	1	1	1
2 ID				
SD				
Group	1	2	2	2

Number of groups

Group(s)	
1	2
1	6
2	7
3	8
4	
5	

Group description:

Here you can configure the experiments and save different information about the animals and the experiment.

You can put in a description of the experiment as well as a description for each animal.

You can also combine animals to a group so the data are summarized later on automatically.



Program Tour - Configuration

FST Analysis

Technical preferences

Experiment Technical preferences

60 208 88 240

Animal 1	Animal 2	Animal 3	Animal 4
Struggle	Struggle	Struggle	Struggle
Swim	Swim	Swim	Swim
Float	Float	Float	Float

Animal 5	Animal 6	Animal 7	Animal 8
Struggle	Struggle	Struggle	Struggle
Swim	Swim	Swim	Swim
Float	Float	Float	Float

Sensitivity 7

Source ☐ Camera ☒ Video file

Video file ☐ Save file

Floating/Active Strugg. - Up area Strugg. - Middle Strugg. - Down

Area 18

Activity 21

Mice		1	2	3
1	Area	3,375	7,5	2,6875
	Activity	17	32	14
2	Area	26,6875	11,1875	21,9375
	Activity	67	38	58

Area

Activity

In this section the filter settings for the detection of the three different behaviors are done.

This only needs to be done once, when the whole system is set up. After that, no filter adjustment is necessary any more.


During the filter configuration the video input and the detected behaviors are shown in real time so it is very easy to adapt the filters until the results meet your expectations.



Program Tour - Acquisition

FST Analysis

Data acquisition during the experiment



In the upper left corner, the video input of the 8 animals is shown. Underneath the video image, the currently detected behaviors are presented in real time for each single animal.

On the right side the summarized data are shown in a table and in diagrams, again for each single animal and also for the animal groups.


The table and the diagrams are updated in dependence of the analysis interval selected in the configuration (in this example every 60 seconds).



Program Tour - Acquisition

FST Analysis

Finished experiment



The table and the diagrams show the data of the entire experiment (360 seconds). In the table there is now also a row that shows the mean values for each behavior.

All data and diagrams can be exported in various file formats.



Program Tour - Data Analysis

FST Analysis

The integrated data analysis offers extensive possibilities to analyze and visualize the data. You can for example group the data of different animals of several experiments and perform intra- and inter-group calculations and comparisons.

This is a list of the provided values and calculations:

Individual analysis

- Min and max values
- Average
- Trend
- Exponential trend
- Polynoms
- p-values of the functions
- Standard deviation
- Standard error
- Upper approximate value

Group analysis

- Min and max values
- Average
- Trend
- Exponential trend
- Polynoms
- p-values of the functions
- Standard deviation
- Standard error
- Upper approximate value
- Standard dev. for single values
- Standard err. for single values
- St. dev. for single animals compared to group average
- St. err. For single animals compared to group average

Group to group analysis

- Min and max values and statistical relevant differences
- Average and statistical relevant differences
- Trend
- Exponential trend
- Polynoms
- Standard deviation for single group values
- Standard error for single group values
- Standard deviation for groups against each other
- Standard error for groups against each other
- Combination of single groups statistics and the inter-group statistic for the calculation of the validity of the inter-broug statistic




Program Tour - Data Analysis

FST Analysis

The following figures show an example for some data of three different groups.

Group list


Use	No	Short description
<input checked="" type="checkbox"/>	1	Dopamin Group
<input checked="" type="checkbox"/>	2	Haloperidol Group
<input checked="" type="checkbox"/>	3	Control Group



This figure shows the exponential trend for floating of animal 1 in group 1



Inter group relations



This figure shows the behaviors for the three different groups and you can easily compare the data. All data and graphs can be exported in different file formats.




Program Tour – Report generator

FST Analysis

The integrated **report generator** makes it possible to create customized and standardized experiment reports in no time.

The reports are fully customizable. You can define what data are presented in which way.

This example report shows the basic data and diagrams of three animals.s






Program Tour – Report generator

FST Analysis

This example report shows more detailed data of one animal.





Add-on products

FST Analysis



Custom made infrared illumination devices improve the results of the video analysis. These units also contain a camera. Each unit can be used with two animals.



An A/D video converter can be used to feed analog video material (either live or taped) into the system.



We offer all inclusive "out-of-the-box" systems that include all devices (computer, camera, cables, illumination etc.) that are necessary to set up an automated FST.



Digital video cameras (FireWire/iLink) are used to feed the video signal into the computer with the FST Analysis software.



A remote control to start and stop the program without the keyboard can be useful in behavioral experiments.




Thank you!

FST Analysis

Thank you very much for your interest in our solution.

We hope we could provide you an useful overview over the features and capabilities of our high throughput Forced Swim Test Analysis Solution „FST Analysis“.

If you have further questions, please do not hesitate to contact us.



Europe:

BIOBSERVE GmbH
Koenigswinterer Str. 418
53227 Bonn
Tel.: +49 228 20 160 20
Fax: +49 228 20 160 29

North America:

BIOBSERVE
2125 Center Ave., Suite 500
Fort Lee, New Jersey 07024
Tel.: 201-302-6083
Fax: 201-302-6062

www.biobserve.com info@biobserve.com