

## Home Cage Activity Monitoring

### *MLog*

#### **Europe:**

BIOBSERVE GmbH  
Koenigswinterer Str. 418  
53227 Bonn  
Tel.: +49 228 20 160 20  
Fax: +49 228 20 160 29

[www.biobserve.com](http://www.biobserve.com)   [info@biobserve.com](mailto:info@biobserve.com)

#### **North America:**

BIOBSERVE  
2125 Center Ave., Suite 500  
Fort Lee, New Jersey 07024  
Tel.: 201-302-6083  
Fax: 201-302-6062

This document and its content is confidential and was created for the solely use of the intended audience. It is not complete unless supported by the underlying analyses and oral presentation of BIOBSERVE GmbH. It may not be reproduced, disclosed or passed on to third parties except with the explicit prior consent of BIOBSERVE GmbH.



## Features I

Home Cage Activity Monitoring - *MLog*

### Measuring activity

The infrared sensors of the MLog system measure the activity of an animal by recording the body heat of that animal and its spatial movement. The activity is measured in all three dimensions in the range of the sensor. Since infrared radiation cannot pass through standard cage plastics, sensors in neighbouring cages do not interfere with each other. The system can be used with any cage size.




### Ten times per second

Each sensor measures the activity of an animal ten times per second (10 Hz), independently of the total number of sensors in use. In competitive systems the data rate decreases with an increasing number of sensors. Although we work with a high data frequency, MLog can run over 24 hours a day 7 days a week and enables you to collect data for your research studies including for phenotyping or to simply monitor the health and welfare of your animals.




### Individual and pooled data

You can divide the animals into groups (e.g. control and test doses) and the software calculates and presents the pooled data in real time. Thus you can compare the results of the different groups even while the experiment is ongoing.



### Real time data analysis and visualization

At the same time, during the experiment the data are also being analyzed and can be presented in user friendly graphs and diagrams. This feature enables you to evaluate the experiment even while it is running so you can decide if you want to continue or stop the experiment and refine the selected parameters.





## Features II

Home Cage Activity Monitoring - MLog

### Further processing of data and graphics

The data, once collected, can be exported to other applications. Further individual data manipulations can then be conducted. All diagrams and graphics which have been generated during the experiment are saved in vector graphic format (wmf) thereby making them easy to integrate and work with in other programs.

#### Export filters for the data:

\*.txt, \*.doc, \*.csv, \*.xls, \*.xml

#### Export filters for the graphs:

\*.bmp, \*.wmf, \*.emf

### State of the art network technology

The sensors communicate with the computer via the TCP/IP protocol, the standard data communication protocol. Thus the data flow is fast and reliable and you can connect several hundreds of sensors to one computer. The sensors are linked with stable standard (Ethernet CAT-5) network cables whereby the same cable transports power and data simultaneously. This makes for very easy wiring.



### Remote control

The MLog system can be accessed via the Internet or the Intranet so making it possible to control an experiment or to analyze data from virtually anywhere in the world. Thus our scientist clients can be based in one location with their experiments being conducted in another, distant location. And teams of scientists in multiple locations can run shared experiments and all be able to access their data.



### Phenotyping integrated in biocontainment systems

MLog can be integrated into biocontainment housing systems such as the BioZone IVC systems (<http://www.biozoneglobal.com>). The ability to remotely control the experiment is ideal since it reduces the necessity to enter the animal area either to perform experiments or to analyse data.





## System Components

Home Cage Activity Monitoring - MLog

The MLog system consists of different hardware components:

**Activity Unit:** contains the IR-sensor and is mounted on the cage

**Central Unit:** manages the data transfer between several Activity Units and the computer


**Cage Holder:** used to mount the Activity Unit on the cage

### Front of the Activity Unit



Each Activity unit gets a unique number, shown in the display at the front. Thus it is very easy to assign the acquired data to each cage.

### Cage holder



The stable cage holder is fixed on the cage. The Activity Unit is plugged on the holder.

RJ45 port to wire the units with stable Cat.5 cables. Power and data go through the same cable, daisy chaining the units.

### Activity Unit

Serial port to connect and control external hardware (e.g. food supplier)

Additional sensor input





## System Architecture

Home Cage Activity Monitoring - *MLog*

There are several options for setting up the MLog system.

In the **standard installation**, the MLog hardware (Activity Units and Central Unit) is installed in the animal room and the MLog software is installed on one computer that is directly connected with the Central Unit. To control the experiment or analyze the data one has to be present in the lab.

In the **enhanced installation**, the MLog software is running on an Application server, allowing one to control and view the experiment from different computers at different locations with the MLog client software. With this solution, different experiments can be performed at the same time and different individuals can simultaneously use the system.




## Standard Installation

Home Cage Activity Monitoring - MLog

All the MLog hardware is located in the animal room and the computer with the MLog software is connected directly to the Central Unit. One experiment can be performed at a time and only at this computer can the experiment be controlled and analyzed.

### Standard option








# Enhanced Installation

## Enhanced option



The MLog hardware is installed as described previously.

However the Central Unit is connected to the MLog Data Server, a computer, that is located near the animal room. This computer is connected via the Intra- or Internet with the MLog Application Server which manages all the configuration and data files. With the MLog client software, it is now possible to log onto the Application Server from any computer connected to the Intra- or Internet. You can control and analyze your experiments from multiple locations or you can, for example, observe an experiment live during a conference call. Even with a PDA or a cellphone you can control the data acquisition. It is also possible to perform multiple experiments at the same time with this set up.



# Program Tour - Configuration/Activity Unit Management

**Configuration - Unit management**

**Activity units**

Auto detect

Use	Unit	Address
<input checked="" type="checkbox"/>	Unit 1	0000
<input checked="" type="checkbox"/>	Unit 2	0000
<input checked="" type="checkbox"/>	Unit 3	0000
<input checked="" type="checkbox"/>	Unit 4	0000
<input checked="" type="checkbox"/>	Unit 5	0000
<input checked="" type="checkbox"/>	Unit 6	0000
<input checked="" type="checkbox"/>	Unit 7	0000
<input checked="" type="checkbox"/>	Unit 8	nnnn

**Groups**

Group 1: Unit 1

Group 2: Unit 2

Description Group1: Control

Description Group2: Hyper

**Activity unit description**

Address: 0000

Type of sensor: Infrared

Description: Control

**Central units**

Number: 0

Current status: no information

Boot option: no information

Central unit Activity units

**Application Server**

☒ Use own IP Address ☐ Auto log in

Address: 192.168.200.105

## Activity units

All activity units in use are shown.

## Activity unit description

You can enter a description for each single activity unit (i.e. for each animal). This may be an animal number or test substance for example.

## Groups

You can group the units and enter a description for each group. Even during the experiment, the group data will be calculated and presented in real time.

## Application Server & Central units

Technical preferences for internal data communication.





# Program Tour - Configuration/Experiment Description

Home Cage Activity Monitoring - MLog

**Configuration - Experiment description**

Experiments	Filename
Experiment: 1 of 1 + - Status: Stopped	demo.mlg
Description	Time and date
Experimenter: <input type="text"/> No.: 0 Description: <input type="text"/>	<input checked="" type="radio"/> real time <input type="radio"/> absolute time <input type="checkbox"/> Use automatic start time

## Experiments

Number of experiments in the system.

## Description

Description of the experiment.

## Filename

Filename of the experiment.

## Time and date

You can select whether real time or absolute time is used for the experiment and you can use the automatic start and stop option.



# Program Tour - Configuration/Experiment Description

**Configuration - Data Analysis**

Activity classes			
Class	from	to	
1	0	2	
2	3	4	
3	5	6	
4	7	8	
5	9	10	
6	11	12	
7	13	14	
8	15	16	
9	17	18	
10	19	20	

**Interval data**

Interval length [hh:mm]: 00:10

☐ Save raw data

**Actogram**

Interval length [hh:mm]: 00:50

Threshold: 1

☒ averaged ☐ percentage: 50

## Activity classes

Division of the ten activity classes.

## Interval data

Time interval of the summary data.

## Actogram


Time interval and threshold of the actogram (also can be changed later).




## Program Tour - Data Acquisition

Home Cage Activity Monitoring - MLog

During the experiment, the data are calculated and presented in real time. Thus it is very easy to evaluate an experiment even while it is ongoing.



For each animal, the current activity level is presented.



Also the activity classes are shown for each unit.

Data							
Time	Date	Activity	Class 1	Class 2	Class 3	Class 4	Class 5
23:00:23	05.08.2003	1	522	18	0	0	
23:01:23	05.08.2003	1	517	20	0	0	
23:02:23	05.08.2003	1	519	31	0	0	
23:03:23	05.08.2003	1	517	19	0	0	
23:04:23	05.08.2003	1	515	34	0	0	
23:05:23	05.08.2003	1	490	43	0	0	
23:06:23	05.08.2003	1	483	52	0	0	
23:07:23	05.08.2003	2	442	73	21	9	
23:08:23	05.08.2003	14	178	119	77	66	3
23:09:23	05.08.2003	3	355	153	33	9	
23:10:23	05.08.2003	3	312	212	14	0	
23:11:23	05.08.2003	4	367	122	32	12	
23:12:23	05.08.2003	1	501	28	6	0	
23:13:23	05.08.2003	2	449	65	13	7	
23:14:23	05.08.2003	1	520	30	0	0	
23:15:23	05.08.2003	1	526	7	0	0	
23:16:23	05.08.2003	1	505	45	0	0	

Of course all data are also presented.




## Program Tour - Data Analysis I

Home Cage Activity Monitoring - MLog


The integrated data analysis offers a wide range of options for analysis and visualisation of the data. You can, for example, group the data of different animals from several experiments and then perform intra- and inter-group calculations and comparisons.

### Graphs



This figure shows the activity rhythm of a rat over four days. You can clearly see the circadian rhythm. MLog allows you to compare the activity data of animals or groups very simply. Just select the units you want to compare in the pop up window and the diagram is redrawn with the selected data. Thus it is very easy to compare different animals or groups. Of course the diagrams can also be printed and saved.

This figure is animated in the Power Point file.



This bar chart shows the time the selected animal spends in each of the ten possible activity classes. The activity classes are a means to divide the activity values into ten groups to provide a more detailed view of the activity seen.




## Program Tour - Data Analysis II

Home Cage Activity Monitoring - MLog

### The actogram

An actogram shows the activity/inactivity of an animal based upon a definable threshold. White symbolizes activity, black inactivity



The actogram is calculated based on two parameters; the activity threshold and the time interval. Both parameters can be changed and the actogram is recalculated in real time. The diagram can be copied and saved in different file formats.

This figure is animated in the Power Point file.




# Thank you!

Home Cage Activity Monitoring - MLog

Thank you very much for your interest in our solution.

We hope we could provide you an useful overview over the features and capabilities of our Home Cage Activity Monitoring Solution „MLog“.

If you have further questions, please do not hesitate to contact us.



**Europe:**

BIOBSERVE GmbH  
Koenigswinterer Str. 418  
53227 Bonn  
Tel.: +49 228 20 160 20  
Fax: +49 228 20 160 29

**North America:**

BIOBSERVE  
2125 Center Ave., Suite 500  
Fort Lee, New Jersey 07024  
Tel.: 201-302-6083  
Fax: 201-302-6062

[www.biobserve.com](http://www.biobserve.com) [info@biobserve.com](mailto:info@biobserve.com)